
Name Scientific Processes - Tools and Measurements


Part A - Count Your Drops

 Take a guess – how many drops of water in one milliliter of water 	

2.	Tools	vou	will	need







*Note, to measure volume with a graduated cylinder, you should take the reading at the lowest point in the curve.

What amount of liquid is in the graduated cylinder pictured? _____How much liquid can your graduated cylinder measure? _____ How much liquid can your beaker measure?

3. To determine the number of drops it takes to make a milliliter, fill you cylinder to 10 ml of water. Carefully add drops using the pipette until you reach 11 ml. Repeat this process 3 times in order to calculate an average.

Trial 1	Trial 2	Trial 3	Average
			-

Part B - The Volume of Solid Objects

1. Solid objects have a volume also (basically the amount of space the object takes up). Volume can be measured in two ways.

For symmetrical objects, volume is simply LENGTH x WIDTH x HEIGHT Use a metric ruler to measure the box below and determine its volume (measure in cm)



The volume of the box