

## Photosynthesis and Cellular Respiration Rates in Plants Using Vernier Probes

**Time:** 55 minutes

**Grade level:** 9-12

### Objectives:

- Students will be able to use CO<sub>2</sub> Gas Sensors to measure the amount of carbon dioxide consumed or produced by a plant during respiration and photosynthesis.
- Students will be able to determine the rate of respiration and photosynthesis of a plant.

### Idaho Achievement Standards:

- 649.01c Use technology and mathematics to improve investigations (data collection and analysis) and communication.
- 652.02a Know that atoms and molecules cycle among the living and nonliving components of the biosphere.
- 652.02b Trace energy flows through ecosystems in one direction, from photosynthetic organisms to herbivores to carnivores and decomposers.
- 653.01b Know that living systems require a continuous input of energy to maintain their chemical and physical organization.
- 653.01c Know that the energy for life is primarily derived from the sun through photosynthesis.
- 653.01d Understand cellular respiration and the synthesis of macromolecules.
- 653.01e Know that chemical bonds of food molecules contain food energy, which is released when the bonds are broken.
- 653.01h Trace how matter cycles and energy flows through different levels of organization of living systems – cells, organs, organisms, communities – and between living systems and the physical environment.

### Materials:

- |  |                                       |                            |
|--|---------------------------------------|----------------------------|
| ▪ Living leaf samples from different types of plants | ▪ Vernier CO <sub>2</sub> Gas Sensors | ▪ Forceps                  |
| ▪ Computers loaded with Logger Pro software          | ▪ 250-mL conical flasks               | ▪ Tape                     |
| ▪ Vernier computer interface                         | ▪ Lamps                               | ▪ Scissors                 |
|  | ▪ Colored cellophane                  | ▪ Paper towels             |
|  | ▪ Ice bath                            | ▪ Colored chalk            |
|  | ▪ Aluminum foil                       | ▪ 500 mL beaker with water |

### Procedure:

- I will engage the students through quickly explaining that they will be conducting experiments on plant leaves to determine the rate of consumption or production of CO<sub>2</sub> using gas sensors. (Time 3 min.)
- I will have the students get into groups of 2-3, come to the materials table and pick up the “Procedures Paper”, “Data Collection/Answer Sheet”, their “Take Home Assignment”, leaf sample (of choice), and 250 mL flask. and let them choose the plant material they wish to work with. (Time 2-3 min.)