

Name: _____

Date: _____

Parts of the Microscope Matching Activity

How do you stop a stinky fish from smelling? Hold its nose.

Clues are listed below. Print the word that matches the clue on the blank line by the clue.



1. _____ Located on the side of the frame, used to adjust the focus of the microscope
2. _____ Supports upper part of the microscope
3. _____ Located on the stage, adjusts the amount of light passing into the slide
4. _____ The bottom of the microscope, used for support
5. _____ Light or mirror that projects light through the diaphragm
6. _____ Connects the eyepiece to the objective lenses
7. _____ Holds the slide in place on the stage
8. _____ The lens or system of lenses in a microscope that is nearest the object being viewed
9. _____ Holds two or more objective lenses and can be rotated to change power
10. _____ Combination of lenses at the viewing end of optical instruments
11. _____ Small platform where the specimen is mounted for examination

Tube	Eyepiece	Stage	Eyepiece
Stage clips	Diaphragm	Objectives	Arm
Adjustment knob	Light source	Base	