

Name: \_\_\_\_\_

## Metric Liquid Volume


Directions: Label the diagram with the proper liquid measures.



1 milliliter



5 milliliters is about  
a teaspoon of liquid



1000 milliliters is 1 liter

Directions: Complete the conversions below.

- 1 liter = \_\_\_\_\_ milliliters
- 5000 ml = \_\_\_\_\_ liters
- 4 liters = \_\_\_\_\_ milliliters
- 9,000 ml = \_\_\_\_\_ liters
- 11 liters = \_\_\_\_\_ milliliters

Directions: Compare the measures below.

- 2 liters \_\_\_\_\_ 50 milliliters
- 5000 ml \_\_\_\_\_ 5 liters
- 15 liters \_\_\_\_\_ 30,000 ml
- 1 liter \_\_\_\_\_ 100 ml
- 1500 ml \_\_\_\_\_ 3 liters

Label the liter container using  
milliliter segments. Fill the  
container with red liquid up to  
500 milliliters.

